

Agency ID:

FACILITY/EQUIPMENT PRESURVEY INFORMATION

Facility name:

Person Contacted:

Survey Date:  (MM/DD/YY)      Surveyor ID:

Manufacturer:       Date of Manufacture:  (month/year)

Model Number:       Serial Number:

A. GENERAL INFORMATION

1.a Facility ID number (FDS):       1.b System ID number (CDRH):

2. Room number or location of radiographic system:

3. Is the operator required to stand in a shielded area while making an exposure? (Y/N)

4. Is a radiation warning label on the control console of the X-ray unit? (Y/N)

5. Are the exposure settings (technique factors) visible to the operator before an exposure is taken? (Y/N)

6. If multiple tubes are controlled by a single exposure switch, is there a visible indication of the tube selected at the control console and at the tube head? (Y/N/X)

7. Is there means to indicate when the beam is perpendicular to the plane of the image receptor? (Y/N)

8. Indicate if each gantry lock functions properly (Y/N/X)      a. Vertical       b. Rotational       c. Along table       d. Across table

9. Does the centering detent work properly? (Y/N/X)

10. Do electrical cables appear to be in good condition? (Y/N)

11. Does the equipment appear to be mechanically stable? (Y/N)

12a. Are up-to-date radiographic technique charts in use in this room? (Y/N)

12b. Are back-up technique charts available for manual mode use of phototimed equipment? (Y/N/X)

13. Are tube rating and tube cooling charts posted? (Y/N)

14. Is a heat unit indicator installed on this unit? (Y/N)

15. Are routine quality control tests performed on this unit? (Y/N)

16. Does a QA or maintenance logbook exist for this unit? (B/N)

MINIATION (Set the SOURCE to LIGHT DETECTOR distance at 100 cm)

1.  $\frac{\text{Total}}{\text{Ambient}} = \frac{\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}}{\phantom{00}} \text{ fc}$

2.  $\frac{\text{Total}}{\text{Ambient}} = \frac{\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}}{\phantom{00}} \text{ fc}$

3.  $\frac{\text{Total}}{\text{Ambient}} = \frac{\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}}{\phantom{00}} \text{ fc}$

4.  $\frac{\text{Total}}{\text{Ambient}} = \frac{\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}}{\phantom{00}} \text{ fc}$

INITIAL SET-UP

1.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  kVp      2.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mA      3.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  sec      or      4.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mAs

REPRODUCIBILITY

1a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	1b. -----	2a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	2b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec
3a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	3b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	4a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	4b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec
5a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	5b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	6a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	6b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec
7a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	7b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	8a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	8b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec
9a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	9b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	10a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR	10b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec

MA LINEARITY

1.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mA      2.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mR      3.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mR      4.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mR      5.  $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$  mR

RADIATION OUTPUT LINEARITY

1. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ kVp (actual)	$\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ kVp (indicated)			
mA station selected	timer station selected	mAs station selected	focal spot (L/S)	measured output
2a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mA	b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	c. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mAs	d. $\boxed{\phantom{00}}$	e. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR
3a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mA	b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	c. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mAs	d. $\boxed{\phantom{00}}$	e. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR
4a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mA	b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	c. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mAs	d. $\boxed{\phantom{00}}$	e. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR
5a. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mA	b. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ sec	c. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mAs	d. $\boxed{\phantom{00}}$	e. $\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}\boxed{\phantom{00}}$ mR

RADIATION OUTPUT LINEARITY (continued)

mA station selected	timer station selected	mAs station selected	focal spot (L/S)	measured output
6a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
7a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
8a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
9a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
10a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
11a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
12a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mA	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	c. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mAs	d. <input type="text"/>	e. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR

TIMER LINEARITY

1.    kVp      2.    mA

Indicated Time	Measured Exposure	Indicated Time	Measured Exposure
3a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR	4a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
5a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR	6a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR
7a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR	8a. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> sec	b. <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> mR

BEAM QUALITY

1.    kVp selected         kVp

2a.    mR      using      0 . 0 mm AL

3a.    mR      using      b.   mm AL

4a.    mR      using      b.   mm AL

5a.    mR      using      b.   mm AL

6a.    mR      using      b.   mm AL

7a.    mR      using      b.   mm AL

ACCURACY

1a. Indicated <input type="text"/> <input type="text"/> <input type="text"/> kVp 3a. <input type="text"/> <input type="text"/> <input type="text"/> kVp	Measured b. <input type="text"/> <input type="text"/> <input type="text"/> kV b. <input type="text"/> <input type="text"/> <input type="text"/> kV	Indicated 2a. <input type="text"/> <input type="text"/> <input type="text"/> kVp 4a. <input type="text"/> <input type="text"/> <input type="text"/> kVp	Measured b. <input type="text"/> <input type="text"/> <input type="text"/> kV b. <input type="text"/> <input type="text"/> <input type="text"/> kV
5. Indicate which kV measurement was calculated      A = Average    E = Effective    P = Peak <input type="text"/>			

kV COMPENSATION (SET AT 80 kVp, KEEP mAs CONSTANT)

Selected mA 1a. <input type="text"/> <input type="text"/> <input type="text"/> 3a. <input type="text"/> <input type="text"/> <input type="text"/>	Measured kV b. <input type="text"/> <input type="text"/> <input type="text"/> b. <input type="text"/> <input type="text"/> <input type="text"/>	Focal Spot c. <input type="text"/> c. <input type="text"/>	Selected mA 2a. <input type="text"/> <input type="text"/> <input type="text"/> 4a. <input type="text"/> <input type="text"/> <input type="text"/>	Measured kV b. <input type="text"/> <input type="text"/> <input type="text"/> b. <input type="text"/> <input type="text"/> <input type="text"/>	Focal Spot c. <input type="text"/> c. <input type="text"/>
---	---	--	---	---	--

SID DETERMINATION

1. Indicated Source to Image Distance (SID)	a. <input type="text"/> <input type="text"/> <input type="text"/> in	or	b. <input type="text"/> <input type="text"/> <input type="text"/> cm	
2. Measured distance from focal spot to tabletop	<input type="text"/> <input type="text"/> <input type="text"/> cm			
Measured distance from tabletop to UTIR film plane	<input type="text"/> <input type="text"/> <input type="text"/> cm			

ACTUAL VERSUS INDICATED FIELD SIZE

1. Does the beam limiting device numerically indicate the field size?	(Y/N) <input type="text"/>
2a. Indicated field size ALONG table?	1. <input type="text"/> <input type="text"/> in      or      2. <input type="text"/> <input type="text"/> cm
2b. Indicated field size ACROSS table?	1. <input type="text"/> <input type="text"/> in      or      2. <input type="text"/> <input type="text"/> cm
3a. Measured light field size ALONG table?	<input type="text"/> <input type="text"/> <input type="text"/> cm
b. Measured light field size ACROSS table?	<input type="text"/> <input type="text"/> <input type="text"/> cm

X-RAY / LIGHT FIELD ALIGNMENT

1a. Measured light field size ALONG table?	<input type="text"/> <input type="text"/> <input type="text"/> cm	b. Measured light field size ACROSS table?	<input type="text"/> <input type="text"/> <input type="text"/> cm
2a. Measured X-ray field size ALONG table?	<input type="text"/> <input type="text"/> <input type="text"/> cm	b. Measured X-ray field size ACROSS table?	<input type="text"/> <input type="text"/> <input type="text"/> cm
3a. Total misalignment ALONG table?	<input type="text"/> <input type="text"/> <input type="text"/> cm	b. Total misalignment ACROSS table?	<input type="text"/> <input type="text"/> <input type="text"/> cm

X-RAY FIELD / UTIR CENTERS COMPARISON

1. Measured misalignment between center of X-ray field and center of test film?	<input type="text"/> <input type="text"/> <input type="text"/> cm
---	---

OPERATION

1. Indicate the type of Positive Beam Limitation (PBL) system in use  
1 = Automatically adjusts the X-ray field size  
2 = Prevents the production of X-rays until manual adjustments are made  
3 = System present but disabled or not used  
4 = System not present
2. Can the collimator be adjusted to a field size smaller than the image receptor? (Y/N)
3. Is there an automatic return to PBL mode when the image receptor or the SID is changed? (Y/N)
4. Is X-ray production prevented at an SID where the PBL system was not intended to operate? (Y/N/X)

PBL SIZING

1. Indicated Source to Image receptor Distance (SID) is indicated as  in or  cm
- 2a. Nominal cassette size in UTIR  by  (ALONG table dimension first) b. Units are IN or CM
- 3a. Measured light field size ALONG table?  cm 3b. Measured light field size ACROSS table?  cm
- 4a. Nominal cassette size in UTIR  by  (ALONG table dimension first) b. Units are IN or CM
- 5a. Measured light field size ALONG table?  cm 5b. Measured light field size ACROSS table?  cm

If a wall bucky is present indicate the following information:

6. Measured distance from frontplate of cassette holder to film plane:  cm
7. Measured distance from source to frontplate of cassette holder:  cm
8. Measured misalignment between center of x-ray field and center of test film from wall bucky:  cm

Fill in for either wall bucky (if present) or 2nd SID for table. If using Wall Bucky - ALONG table = VERTICAL

9. Indicated Source to Image receptor Distance (SID) is indicated as  in or  cm
- 10a. Nominal cassette size in Bucky  by  (ALONG table dimension first) b. Units are IN or CM
- 11a. Measured light field size ALONG table?  cm 11b. Measured light field size ACROSS table?  cm
- 12a. Nominal cassette size in Bucky  by  (ALONG table dimension first) b. Units are IN or CM
- 13a. Measured light field size ALONG table?  cm 13b. Measured light field size ACROSS table?  cm

COMMENTS AND OBSERVATIONS